


# Filter Kit for **AeroEclipse\***

Breath Actuated Nebulizers



- ✓ Filters fugitive emissions (aerosolized medication)<sup>1</sup>, viruses<sup>2</sup> and bacteria<sup>2</sup> from the patient's exhaled breath
- ✓ Reduces exposure of healthcare providers to potentially harmful emissions

The filter kit is a single patient, single use accessory to the **AEROECLIPSE\*** Breath Actuated Nebulizer (BAN).

Please contact TMI to find out if the filter kit is available in your market.

Filter Kit shown with **AEROECLIPSE\*** II BAN

### Aerosol Performance:

Aerosol Output<sup>1</sup> of **AEROECLIPSE\*** BAN with Filter Kit compared to BAN alone  $\geq 90 \%$

Aerosolized Medication Emissions<sup>1</sup> when assembled with **AEROECLIPSE\*** BAN  $\leq 1 \%$

### Filtration Performance:

Viral Filtration Efficiency<sup>2†</sup>  $> 99.9970 \%$

Bacterial Filtration Efficiency<sup>2†</sup>  $> 99.9992 \%$

† The bacterial and viral filtration efficiency performance of the filtration media has been tested by the material's manufacturer in accordance to MIL-M-36954C Military Specification, Mask, Surgical, Disposable.

1. TMI, Data on file. 2. Technical Data Sheet, 14 August 2008, on file TMI. T15010-003 Rev B \* trade marks and registered trade marks of Trudell Medical International (TMI). © TMI 2019, 2020. All rights reserved.



**TRUDELL MEDICAL**  
INTERNATIONAL

+1-866-510-0004 Fax: +1-888-314-7233  
customerservice@trudellmed.com **trudellmed.com**

[aeroeclipse.com](http://aeroeclipse.com)



USE ONLY ONCE